

Clinisept+

Cleans • Calms • Cares

Clinisept+ is a next generation antimicrobial skin cleanser for aesthetic and cosmetic treatments.

It offers a completely new approach to professional skin cleansing that maximises protection and optimises patient outcomes





Clinisept+

The information contained herein has been compiled for Total Therapy Training and the British Acupuncture Federation. It contains information. The information provides a factual review of Clinisept+ performance: this exceeds the terminology that is permitted to be used on pack and in advertising to describe the efficacy of Clinisept+, hence the information is not for dissemination beyond the persons to whom it has been sent and should be treated as confidential.

I would be very pleased to host a Zoom meeting to tell you more about Clinisept+ and answer any questions that you may have regarding our technology, its history, method of action and use. There are also many leading practitioners who would be pleased to discuss their experiences of using Clinisept+.



Director
Clinical Health Technologies Ltd
07971 173303

All of the information contained herein is confidential and should not be disseminated beyond the personnel to which it has been sent.

Clinisept+

A gamechanger in skin cleansing - and recovery



Unparalleled protection against complications



Unprecedented skin compatibility



Optimum conditions for rapid skin recovery



Clinisept⁺

Technology: A unique combination of efficacy and safety

EFFICACY

Bactericidal

Fungicidal

Virucidal

Sporicidal

**High log reductions
(99.9999%) in seconds**



SAFETY

Skin neutral pH

Non-irritant

No CHX, alcohol or parabens

Non-cytotoxic

**Causes no irritation or
skin damage**

A high performance antimicrobial solution that is non cytotoxic and leaves no chemical burden on the skin, thereby providing the optimum environment for rapid cell regeneration

Clinisept⁺

In aesthetics

- ***A completely new concept in aesthetic medicine:***
- ***Can be used before, during and after a procedure***
- ***Causes no skin irritation, stinging or discomfort***



Clinisept+

“

Clinisept+ forms a ‘perfect environment’ in which skin regenerates rapidly due to being burdened by neither pathogens or chemistry. It has completely changed the experience of my patients.

”

Paul Baguley

MB.ChB.BSc.LMCC.FRCPC.FRCPCpaeds FRCSC.FRCSCPlas.FRCS(Ed) Plas. FRSM

Clinisept+

Used, and endorsed by UK's leading aesthetic practitioners - including the JCCP

FACE
FACIAL AESTHETIC CONFERENCE AND EXHIBITION

IAPCAM
International Association for Prevention of Complications in Aesthetic Medicine

baon
British Association of Cosmetic Nurses

The ROYAL SOCIETY of MEDICINE

British College of Aesthetic Medicine

IAAGSW
International Association of Aesthetic Gynaecology and Sexual Wellbeing

GREAT BRITISH ACADEMY OF AESTHETIC MEDICINE

BAHRS
BRITISH ASSOCIATION OF HAIR RESTORATION SURGERY

CCR EXPO
CLINICAL + COSMETIC + RECONSTRUCTIVE

Aesthetics
Conference & Exhibition

aesthetic
— MEDICINE LIVE —



"I'm seeing three to four weeks recovery in 10 days"

Mr Paul Baguley
MB ChB BSc LMCC MD FRPC FRCS FRCSed Plas

"I trust and use Clinisept & it delivers incredible results"

Mr Kambiz Golchin
MMed Sci, FRCS(ORL-HNS)



WINNER 2017
Aesthetics Awards
THE BARRY KNAPP AWARD FOR PRODUCT INNOVATION OF THE YEAR,

www.cliniseptplus.com

CLINICAL HEALTH TECHNOLOGIES




vs traditional skin cleansing chemistries: unsurpassed levels of protection and better skin compatibility

	STRENGTH	ANTI-MICROBIAL	SPORICIDAL	NO RESISTANCE/TOLERANCE CONCERNS	NON-IRRITANT TO SKIN, EYES, MUCUS MEMBRANES
2% Chlorhexidine Gluconate	20,000ppm	✓	✗	✗	✗
20% Chlorhexidine Gluconate	200,000ppm	✓	✗	✗	✗
70% Alcohol	700,000ppm	✓	✗	✗	✗
10% Povidone Iodine	100,000ppm	✓	✗	✗	✗
	250ppm	✓	✓	✓	✓

Clinisept+

Additional in use benefits



Reduced patient discomfort



A note about Clinisept+'s efficacy

Clinisept+ is so effective because it utilises an oxidising, rather than toxic, method of action to clean the skin. This enables it to pass all relevant EN standards in seconds.

For example, EN14476 (the recognised measure of efficacy against viruses) requires a 4 log reduction (99.99%) to be achieved within a 5 minute contact time. Clinisept+ achieves a 5 log reduction (99.999%), within 15 seconds. (Test data available).

Unlike chlorhexidine, because of its oxidising method of action, Clinisept+ is also not subject to antimicrobial resistance: IE pathogens cannot become immune to it.



Clinisept+

Cleans • Calms • Cares

No other product combines similar levels of skin hygiene with complete skin and tissue compatibility.

Clinisept+ is used and trusted by the world's leading cosmetic professionals.

“Clinisept has been a great substitute for chlorhexidine for me in both my injectable and laser based procedures. I have noticed reduced skin irritation and erythema with its use versus chlorhexidine both during and after treatments.”

Dr Benji Dhillon

Aesthetic Surgeon. Phi Clinic

“Clinisept is my cleanser of choice for preparing my patients skin in advance of their skin analysis, injectable and laser treatments.”

Miss Sherina Balaratnam

MBBS, MRCS, MSc (UCL)
Medical Director, S-Thetics.